



DIAGNOSTICS

The Spirometry Source

## Astra Product Matrix

	ASTRA 100	ASTRA 200	ASTRA 300
Hand-held Portability	X	X	X
27 Forced Parameters of FVC	-	X	X
NLHEP Parameters	X	X	X
9 Relaxed Parameters of VC	-	X	X
Flow Trigger Tidal Breathing	X	X	X
Bronchodilator Testing	-	X	X
Bronchial Challenge Testing	-	X	X
1500+ Patient Storage	-	-	X
100 Patient Storage	X	X	-
Ethnic Correction %	-	-	-
Ethnic correction (NHANES III)	X	X	X
NLHEP Compliant	-	-	X
5 Button Keyboard	X	X	-
Touchscreen LCD	-	-	X
14 Predicted Sets (Selectable)	-	X	X
4 Interpretations (Selectable)	-	X	X
On-screen Graphics	-	-	X
OSHA Firmware Module	-	-	O
Firmware Personalization	-	X	X
NiMH Battery Kit	O	O	O
USB to Printer	X	X	X
USB to PC	-	X	X
Bluetooth to PC Module	-	O	O
Download test data to PC	-	X	X
Interactive PC Testing	-	X	X
AstraPro Spirometry Software	-	X	X
PIN Protection	X	X	X
Custom Reporting	-	X	X
CAL Check Mode	X	X	X
Flex Pulse Oximetry Module	-	O	O
Upgradeable Firmware	X	X	X

X Standard Feature  
O Optional  
- Feature not included